

Cervical Collar (Hard)



Myofascial Pain



Neck Fracture



Neck Strain



Whiplash Injury



Comfosmart Cervical Collar Hard & Hard Height Adjustable is a premium neck brace designed to provide rigid support and customizable height adjustment for individuals with severe neck injuries or conditions requiring strict neck stabilization. With its durable construction and height adjustability feature, this collar offers maximum stability and optimal comfort. Available in sizes S/M/L/XL, it ensures a precise fit for individuals with varying neck sizes, promoting effective healing and recovery.

Size	Hard	Height Adjustable	Height Adjustable (Premium)
Small	RA5110S	RA5112S	RA5114S
Medium	RA5110M	RA5112M	RA5114M
Large	RA5110L	RA5112L	RA5114L
XL	RA5110XL	RA5112XL	RA5114XL



For Supporting,
Immobilizing or
Adjusting
the Neck in
Required Position

Features:

1. **Rigid Support and Immobilization:** The Cervical Collar Hard (Height Adjustable) provides rigid support and immobilization to the cervical spine, limiting motion and reducing stress on the neck muscles and ligaments for optimal healing.
2. **Height Adjustability:** The collar features a height-adjustable design, allowing for personalized positioning to accommodate individual anatomical variations and comfort preferences.
3. **Durable Construction:** Crafted from high-quality materials, this collar offers exceptional durability and longevity, ensuring reliable support throughout your recovery.
4. **Customizable Fit:** Available in sizes S/M/L/XL, the collar can be adjusted to fit different neck sizes, providing a secure and comfortable fit tailored to individual needs.
5. **Easy Application:** The collar is designed for easy and straightforward application, making it suitable for both healthcare professionals and individuals seeking self-care.

Indications:

- Severe cervical spine injuries, such as fractures or dislocations
- Post-operative immobilization after cervical spine surgery
- Traumatic neck injuries with significant instability
- Cervical muscle strains or sprains requiring rigid support.
- Chronic neck pain with limited mobility and stability

Instructions of Use:

1. Select the appropriate size (S/M/L/XL) based on your neck measurements.
2. Place the collar around your neck, ensuring that the chin rests comfortably on the collar's support.
3. Adjust the height of the collar by sliding it up or down to align with the base of your skull, ensuring a snug fit and proper immobilization.
4. Secure the collar by fastening the Velcro straps at the back of your neck, adjusting the tightness to achieve a secure fit without excessive pressure or discomfort.
5. Verify that the collar provides the necessary support and immobilization required for your condition, limiting neck movement as instructed by your healthcare professional.
6. Wear the collar for the prescribed duration and as directed by your healthcare professional, adhering to their guidance on when to remove the collar for activities such as bathing or sleeping.
7. Regularly check your skin for any signs of irritation or pressure sores. If discomfort, redness, or skin problems occur, consult a healthcare professional for further evaluation and adjustment.

Note: The Cervical Collar Hard is intended for temporary use. Follow the recommendations of your healthcare professional regarding the duration of collar wear and any additional instructions specific to your condition.